STATIM 2000/5000 G4

CASSETTE AUTOCLAVE™

- Operator's Manual

SciCan

| Your Infection Control Specialist™

Contents

1. Introduction	4	6. Storing and Retrieving Cycle Information	32
2. Important Information	5	Retrieving Cycle Information Using the Touchscreen	
Disclaimers		Retrieving Cycle Information Using the USB Data Back up	
Drying Performance		7. Printing Cycle Information	33
Unit Overview - STATIM 2000 G4		Connecting a Printer	
Unit Overview - STATIM 5000 G4		Adjusting your Print Settings	
Touchscreen Overview		External Printers and Specifications	
Set-Up Menu Overview		8. Maintaining your STATIM	34
3. Setting up your STATIM	11	Cleaning the Cassette	
Positioning and Powering your Unit		Cleaning the Water Reservoir Filter	
Levelling your Unit		Cleaning the Water Reservoir	
Connecting the Waste Bottle		Cleaning the Exterior Surfaces	
Filling the Water Reservoir		Changing the STATIM 2000 G4 Air Filters	
Priming the Pump		Changing the Bacteria-Retentive Air Filter	
Setting the Time		Replacing the Cassette Seal	
Setting the Date		Maintaining Fluid Levels	
Setting the Language		Reading Water Quality	
Assigning the Unit Identifier Number		Preventative Maintenance Schedules	
Creating a User ID and PIN		Shipping the Unit / Draining the Reservoir	
Setting Up Process Enforced Usage		9. Troubleshooting your STATIM	40
Changing the Touchscreen Display Themes		10. Test Protocols	44
Adjusting the Screensaver Delay		11. Ordering Spare Parts	46
Adjusting the Screen Contrast		12. Warranty	47
Turning the Button Sound ON or OFF		13. Specifications	48
Adjusting the Button Beep Volume			
Connecting to a Network			
4. Using Cassettes and Preparing Instruments	19		
Using the STATIM 2000 G4 Cassette			
Using the STATIM 5000 G4 Cassette			
Using Drying Plates with the STATIM 5000 G4			
Preparing and Loading Instruments			
Using Biological and Chemical Indicators			
Instrument Weight Guide			
5. Using your STATIM	24		
Selecting a Cycle			
Running a Cycle			
Stopping a Cycle			
Cycle Specifications and Graphs			

STATIM Cassette Autoclave and STATIM are registered trademarks and STAT-DRI, Your Infection Control Specialist, and DriTec logos are trademarks of SciCan Ltd. All other trademarks referred to in this manual are the property of their respective owners.

For all service and repair inquiries:

In Canada 1-800-870-7777
United States: 1-800-572-1211
Germany: +49 (0)7561 98343 - 0
International: (416) 446-4500
Email: techservice.ca@scican.com

Manufactured by:

SciCan

1440 Don Mills Road,
Toronto ON M3B 3P9
CANADA

Phone: (416) 445-1600
Fax: (416) 445-2727
Toll free: 1-800-667-7733

EU Representative

SciCan GmbH
Wangener Strasse 78
88299 Leutkirch
GERMANY

Tel.: +49 (0)7561 98343 - 0
Fax: +49 (0)7561 98343 - 699

SciCan Inc.


701 Technology Drive
Canonsburg, PA 15317
USA

Phone: +1 724 820 1600
Fax: +1 724 820 1479
Toll free: 1-800-572-1211

SciCan Medtech

Alpenstrasse 16
CH-6300 ZUG
SWITZERLAND
Phone: +41 (0) 41 727 7027
Fax: +41 (0) 41 727 7029

1 Introduction

Congratulations on selecting the STATIM® Cassette Autoclave. We are confident that you have purchased the finest equipment of its type. The STATIM is a compact, countertop unit that features a number of sterilizing cycles designed to meet your needs and suitability for steam sterilization. This STATIM G4 Cassette Autoclave is fully compliant with EN13060.

The details of installing, operating and maintaining your STATIM are all contained within this operator's manual. To ensure years of safe, trouble-free service, please read these instructions before operating this unit and keep them for future reference. Operational, maintenance and replacement instructions should be followed for the product to perform as designed. Contents of this manual are subject to change without notice to reflect changes and improvements to the STATIM product.

The STATIM is suitable for the sterilization of dental and medical instruments designed to withstand steam sterilization. The STATIM has not been designed to sterilize liquids, cloth loads, bio-medical waste or materials not compatible with steam sterilization. The processing of such loads may result in incomplete sterilization and / or damage to the autoclave. For more information about instrument suitability for steam sterilization, consult the instrument manufacturers' reprocessing instructions.

2 Important Information

2.1 Disclaimers

Use only steam-process distilled water in your STATIM. Deionized, demineralized, or specially filtered water should not be used. Never use tap water.

Do not permit any person other than certified personnel to supply parts, service or maintain your STATIM. SciCan shall not be liable for incidental, special or consequential damages caused by any maintenance or services performed on the STATIM by a third party, or for the use of equipment or parts manufactured by a third party, including lost profits, any commercial loss, economic loss, or loss arising from personal injury.

Never remove the cover of the unit and never insert objects through holes or openings in the cabinetry. Doing so may damage the unit and / or pose a hazard to the operator.

IMPORTANT: Follow local guidelines governing verification of the sterilization procedure.

Drying Performance

The STATIM 2000 G4 and 5000 G4 have been designed to provide a complete sterilization solution for your unwrapped and wrapped instrument needs: rapid sterilization balanced with rapid drying, by using the SciCan Dri-Tec Drying Technology.

The STATIM 2000 G4 uses convection heat to dry instruments by utilizing the remaining heat in the system after the sterilization phase. Heat is captured and released in the cassette to rapidly dry a properly loaded STATIM cassette.

The STATIM 5000 G4 uses heat generated by the sterilization phase which is absorbed by the drying plates. The heat is transferred from drying plates directly to the load resulting in accelerated, rapid drying of a properly loaded STATIM cassette.

Please refer to this operator's manual for instructions on the proper arrangement of instruments in the cassette and the use of drying plates (STATIM 5000 G4). By carefully following these directions on how to properly load the cassette chamber, rapid drying of the load will be achieved.

2 Important Information 2000 G4

2.2 STATIM 2000 G4 — Unit Overview

- 1** LCD touchscreen
- 2** reservoir cap/
water filter
- 3** USB port
- 4** power switch
- 5** power cable port
- 6** leveler feet
- 7** Ethernet port
- 8** exhaust tube port
- 9** compressor
- 10** cassette
- 11** biological filter



Figure 1

- 12** RS232 port

The following symbols appear in the margins of this book.

A potential hazard to the operator.

A situation which may lead to a mechanical failure.

Important information

The following symbols appear on the unit:

Caution: Hot Surface and/or Hot Steam

Caution: Refer to manual for details








Caution: Risk of electrical shock. Disconnect supply before servicing

Steam distilled water only

When you receive your STATIM 2000 G4, the items listed below will be included. If any of the items are missing, contact your dealer immediately.

	Cassette Tray and Lid
	Instrument Rack
	Waste Bottle
	Bottle Lid Fitting
	Tube Mounting Hardware
	Power Cord

	Operator's Manual
	Exhaust Tube
	STAT-DRI
	P.C.D. + 20 SciCan Sterilization Emulators
	USB Memory Stick

2 Important Information 5000 G4

2.2 STATIM 5000 G4 – Unit Overview

- 1** LCD touchscreen
- 2** reservoir cap/
water filter
- 3** USB port
- 4** power switch
- 5** power cable port
- 6** leveler feet
- 7** exhaust tube port
- 8** biological filter
- 9** cassette
- 10** RS232 port
- 11** Ethernet port


Figure 2

The following symbols appear in the margins of this book.

A potential hazard to the operator.

A situation which may lead to a mechanical failure.

Important information

The following symbols appear on the unit:

Caution: Hot Surface and/or Hot Steam

Caution: Refer to manual for details







Caution: Risk of electrical shock. Disconnect supply before servicing

Steam distilled water only

When you receive your STATIM 5000 G4, the items listed below will be included. If any of the items are missing, contact your dealer immediately.

	Cassette Tray and Lid
	Unwrapped Instrument Rack
	Waste Bottle
	Bottle Lid Fitting
	Tube Mounting Hardware
	Rack with Drying Plates

	Power Cord
	Operator's Manual
	Exhaust Tube
	STAT-DRI
	P.C.D. + 20 SciCan Sterilization Emulators
	USB Memory Stick

2.4 Touchscreen Overview




Figure 3

2.5 Set Up Menu Overview

2.5 Set Up Menu Overview con't

3 Setting up your STATIM

3.1 Positioning and Powering your Unit

Positioning your Unit

There are several factors that may affect the performance of your STATIM. Please review these factors, and select a suitable location in which to install the unit.

- **Temperature and Humidity**

Avoid installing your STATIM in direct sunlight or close to a heat source (e.g. vents or radiators). The recommended operating temperatures are 15-25°C with humidity of 25-70%.

- **Spacing**

The vents and openings on your STATIM should remain uncovered and unobstructed. Leave a minimum of 50 mm between the top, sides and back of the unit and any wall or partition. For more detailed information on clearances, see Specifications.

- **Venting**

Your STATIM should be operated in a clean, dust free environment.

- **Work Surface**

Your STATIM should be placed on a flat, level, water-resistant surface. Never install and operate the unit on a sloped surface.

- **Electromagnetic Environment**

Your STATIM has been tested and meets applicable standards for electromagnetic emissions. While your unit does not emit any radiation, it may itself be affected by other equipment which does. We recommend that your unit be kept away from potential sources of interference.

- **Electrical Connection**

To power your unit, use properly grounded and fused power sources with the same voltage rating as indicated on the label at the back of your STATIM. Avoid multiple outlet receptacles. If using a surge suppressor power bar, connect only one STATIM to it.

Powering your Unit

To power your STATIM, connect the power cord to the A.C. inlet receptacle at the back of the unit. Ensure the power switch is in the OFF position and connect the unit to your power source.

3 Setting up your STATIM


3.2 Levelling your Unit

When placing your unit on a counter top, ensure it is stable and that all four feet are securely in contact with the counter surface. This will prevent the unit from moving freely. Next, use the level indicator bubble function in the settings menu to adjust the three leveler feet so that the unit drains properly. To access the level indicator bubble from the main screen, follow these steps:

2. Scroll to and select.

3. Adjust the leveler feet to move the bubble. Position the bubble in the bottom right quadrant of the target. This will ensure that the unit drains properly. Press STOP to exit and return to the cycle select menu. When the unit is properly leveled, the bubble level will change from red to green.

3 Setting up your STATIM

3.3 Connecting the Waste Bottle


The waste bottle **1** is used to collect the wastewater after it has been converted to steam and then drained from the cassette. To connect the waste bottle to the STATIM, follow these steps (see figure 4):

1. Insert the exhaust **2** tube into the fitting **3** on the back of the unit and pull gently to confirm a tight fit.
2. Cut the tube to length and slide the waste bottle fitting **4** into place.
3. Place the free end of the tube into the hole in the lid of the waste bottle and hand-tighten the fitting. Do not coil the exhaust tube.
4. Unscrew the lid and copper condenser coil **5** assembly from the waste bottle. The lid and coil should come out together.
5. Fill the waste bottle with water to the MIN line and replace the lid and copper condenser assembly. Empty the waste bottle often to avoid unpleasant odors and discoloration of the contents. (A low-level disinfectant, prepared according to the manufacturer's instructions, may be added to the waste bottle to remedy this situation). As a minimum, empty the waste bottle each time you refill the reservoir.

6. Place the waste bottle near the unit. Store the bottle below the unit. The tube can be routed through a hole, (8 mm in diameter) in the counter-top and secured with the provided nylon clamps.

3 Setting up your STATIM

3.4 Filling the Water Reservoir

When filling the reservoir, ensure you only use steam processed distilled water containing less than 5 ppm total dissolved solids (having conductivity of less than 10 $\mu\text{S} / \text{cm}$).

The impurities and additives in other water sources will cause an error reading on the LCD. To fill the reservoir, follow these steps (see Figure 5):

1. Remove the reservoir cap **1**.
2. Pour steam-process distilled water into the reservoir until almost full (a maximum of 4L). Use a funnel to avoid spillage.
3. Replace and secure the cap.


Figure 5

3.5 Priming the Pump

To prime the STATIM pump, follow these steps:

1. Move the unit to the edge of the work surface. The front leveler feet should be approximately 12 mm from the edge.
2. Lift the front left corner of the unit upward and remove the drain tube **2** from the clip located on the underside of the unit.
3. Pull the drain tube outward so the free end can be positioned over a water container.
4. Fill the reservoir with steam-process distilled water.
5. Remove the plug **3** from the end of the drain tube and allow water to drain from the tube into a container for 30 seconds. When the water flows in a steady stream, replace the stopper.
6. Lift the front left corner of the unit upward and reinsert the tube into the clip on the underside of the unit. Push the excess length of tubing back into the space provided.


Figure 6

Make sure the plug in the drain tube is secure.

3 Setting up your STATIM

3.6 Setting the Time

1. 
2. Scroll to and select.
3. From the TIME screen, use the keypad to set the time. Press to save and  to return to the Setup menu.
4. To change your unit to display 12-hour time format (24-hour time format is the default setting), go to the Setup menu and use   to scroll to TIME 12/24, select it and toggle to 12. Press  to save and return to the Setup menu.
5. To activate daylight savings time (DST), which is recommended if you have connected your unit to a network, go to the Set up menu and use   to scroll to DST ON/OFF and select. Use   to toggle DST ON or OFF and press the  to save and return to the Setup menu.

3.7 Setting the Date

1. 
2. Scroll to and select.
3. From the DATE screen, use the keypad to set the date. Press to save and  to return to the Setup menu.
4. To change the format in which the date appears, return to the Setup menu and use   to scroll to DATE FORMAT. Select it, and follow the prompts to have the date displayed in the desired format. Press  to save and return to the Setup menu.

3.8 Setting the Language

The messages displayed by your STATIM can be presented in a number of different languages. To change the current language, follow these steps:

1. 
2. Scroll to and select.
3. From the LANGUAGE screen, press   to scroll through the list of languages. When you have found the desired language, press  to save your selection and return to the Setup menu.


3 Setting up your STATIM

3.9 Assigning Unit Identifier Number

1. 
2. Scroll to and select.
3. Using the keypad, select a maximum of 3 digits to be used as the unit's identifier number. Press to save and  to return to the Setup menu.

3.10 Creating a User ID and PIN

1. 
2. Scroll to and select.
3. From the SETUP PIN screen, you can assign up to four PINs. Select one of the User icons to assign a PIN.
4. Using the keypad, assign a PIN of up to four digits and press to save and  to move to the confirmation screen.

5. If all of the information presented in the confirmation screen is correct, press OK to be returned to the PIN USER screen. To make a correction, select the PIN User you want to change and repeat the process described above.

3.11 Setting Up Process Enforced Usage

When process enforced usage is activated, users are required to enter a PIN both at the beginning and at the end of a cycle. For Process Enforced usage to function, User IDs and PINs must first be assigned. To set up User ID and PINs, refer to section 3.10 Creating a User ID and PIN. To activate process enforced usage, follow these steps:


1. 
2. Scroll to and select.
3. Use   to toggle Process Enforced function ON or OFF. Press  to save your selection and return to the Setup menu.

NOTE: Any user can stop a cycle and remove the cassette even with Process Enforced usage ON. However, the cycle data will record that an unauthorized user has stopped the cycle and/or removed the cassette.

3 Setting up your STATIM

3.12 Changing the Touchscreen Display Themes

Your STATIM G4 touchscreen themes (i.e., icons and background colours) can be changed to one of the preset options or you can upload additional themes, as they become available from SciCan, using the USB port. To change themes follow these steps:

1. 
 2. Scroll to **Themes** and select.
 3. From here, you can either select **Change Theme** for a menu of preloaded themes or **Upgrade Theme** to access a new theme that can be loaded using the USB port.
 4. In the **Change Theme** screen, use   to scroll through your available options. As you scroll, each theme will display on the touchscreen. Press  to select your theme and return to the Setup menu.
 5. To upgrade a theme available from SciCan, download the theme onto your computer's desktop and save the files onto a portable USB storage device. Insert the device into your STATIM's USB port and, from the UPGRADE THEME screen, press 
- 5.1** The unit will load the files directly from the USB storage device. Do not remove the USB storage device while files are loading (this could take as long as 10 minutes). When it is complete, the screen will display the 'Done' message. This new theme will now be accessible from your THEMES menu.

5.2 Press  to select this theme and return to the Setup screen.

3.13 Adjusting the Screensaver Delay

To change the length of time before your inactive touchscreen activates the screensaver, follow these steps:

1. 
2. Scroll to **Screensaver** and select.
3. Use   to scroll through your time options. When you have found the amount of time you require, press it. Press  to save and return to the Setup menu.

3 Setting up your STATIM

3.14 Adjusting the Screen Contrast

The STATIM G4 touchscreens are calibrated for the lighting condition of most sterilization centres. Should you need to adjust the contrast for your office, follow these steps:

1. 
2. Scroll to **LCD Contrast** and select.
3. Use   to scroll through your contrast options. When you have found the contrast you require, press it. Press  to save and return to the Setup menu.

3.15 Turning the Button Sound ON or OFF

The STATIM G4 is preset to beep when a button is pressed. If you would like to turn the button sound off, follow these steps:

NOTE: Turning OFF the button sound does NOT turn off other alarms and cycle notification beeps.

1. 
2. Scroll to **Beep ON/OFF** and select.
3. Use   to scroll through your ON or OFF options and select it by pressing it. Press  to save and move back to the Setup menu.

3.16 Adjusting the Button Beep Volume

If you would like to adjust the beep volume, follow these steps:

1. 
2. Scroll to **Beep Volume** and select.
3. Use   to scroll through to scroll through the volume settings. Select the one you want by pressing it. Press  to save and move back to the Setup menu.

3.17 Connecting to a Network

The STATIM G4 has a 10/100Base-T Ethernet port located at the back of the unit. To connect your STATIM to a network consult the separate manual titled STATIM 2000/5000 G4 – Setting Up and Using Your Web Portal.

4 Using Cassettes and Preparing Instruments


4.1 Using the STATIM 2000 G4 Cassette

When removing the cassette after a cycle, exercise caution as the metal areas will be hot and the cassette may contain hot steam.

Opening the Cassette:

1. Hold the cassette handle with your thumbs facing inward on the cassette latch.
2. Push downward on the cassette latch.
3. Raise the cassette lid upwards and disengage the hinge.
4. Rest the lid on its outer surface.

Closing the Cassette:

1. Align the hinge tab on the cassette lid with the hinge slot on the rear of the bottom tray.
2. As you begin to close the lid, the hinge tab and slot will engage.

Inserting the Cassette into the STATIM G4 2000:

1. Place the end of the cassette into the unit.
2. Gently push inward until you hear a “click” sound and check the LCD icon for change from  to .

Never force the cassette into the STATIM as the interior components could be damaged.

NOTE: The main menu screen will display  if the cassette is not properly inserted in the unit.

Removing the Cassette:

1. Grasp the handle with two hands and pull away from the unit.
2. Pull the cassette clear of the unit and set down on a firm surface.

Disengaging the Cassette:

When not in use, the cassette should be disengaged. To disengage the cassette, grasp the handle and pull the cassette out until there is a 15 mm to 20 mm ($1/2$ to $3/4$ ") gap between the front of the STATIM and the cassette handle.

Using the STAT-DRI Drying Agent

Treatment of the interior surfaces of the cassette with the STAT-DRI drying agent, provided with your unit, will enhance the drying process. (Replacement bottles are available from SciCan, order number 2OZPLUS, 8OZPLUST, 32OZPLUS).

4 Using Cassettes and Preparing Instruments

4.2 Using the STATIM 5000 G4 Cassette

When removing the cassette after a cycle, exercise caution as the metal areas will be hot and the cassette may contain hot steam.

Opening the Cassette:

1. Push the carry handle **2** into the open position.
2. Put your hands on either side of the cassette handle **1**.
3. Insert your forefingers in the slots and place your thumbs on the thumb pads.
4. Press down with your thumbs and pull up with your forefingers until the lid opens.
5. Raise the cassette lid and disengage from the tray. Rest the lid on its outer surface.

Closing the Cassette:

1. Align the hinge tab on the lid with the hinge slot on the tray.
2. As you begin to close the lid, the hinge tab and slot will engage.
3. Place the carry handle **2** into the closed position.


Inserting the Cassette into the STATIM 5000 G4:

1. Hold the cassette handle in one hand and the carry handle in the other as shown in Figure 8.
2. Place the end of the cassette into the unit and drop the carry handle into its closed position.
3. Gently push inward until you hear a “click” sound and check the LCD icon for change from  to .

Never force the cassette into the STATIM as the interior components could be damaged.

NOTE: The main menu screen will display  if the cassette is not properly inserted in the unit.

4 Using Cassettes and Preparing Instruments

Removing the cassette:

1. Grasp the cassette handle with one hand and pull out from the unit.
2. As the cassette emerges from the unit, grasp the carry handle with your free hand and lift it upwards.
3. Pull the cassette clear from the unit and set down on a firm surface.

Disengaging the cassette

When not in use, the cassette should be disengaged. To disengage the cassette, grasp the handle and pull the cassette out until there is a 15 mm to 20 mm ($1/2$ to $3/4$ ") gap between the front of the STATIM and the cassette handle.

4.3 Using Drying Plates with the STATIM 5000 G4

The STATIM 5000 G4 cassette comes with two racks, one with drying plates affixed and another without. Use the rack with drying plates to enhance the drying process for wrapped instruments.


Figure 9

4.4 Preparing and Loading Instruments

Before loading any instruments into the STATIM, consult the manufacturer's reprocessing instructions.


Clean Instruments

Clean and rinse all instruments before loading them into the cassette. Disinfectant residues and solid debris may inhibit sterilization and damage the instruments, the cassette, and the STATIM. Lubricated instruments must be wiped thoroughly and any excess lubricant should be removed before loading.

Unwrapped Instruments

Arrange unwrapped instruments on the rack spreading them out as evenly as possible.


4 Using Cassettes and Preparing Instruments

Wrapped Instruments (single wraps)

Place the instruments into single layer autoclave bags according to the manufacturer's instructions. Orient the instrument rack in the cassette to ensure that wrapped instruments rest approximately 6 mm / 0.25" above the cassette base. Place the wrapped instruments on the rack and arrange them to avoid overlap. Ensure that all wrapped loads are dry before handling and /or storage to maintain sterility.


For STATIM 2000 G4:

The use of cloth wraps in the STATIM is not recommended.

SciCan recommends the use of paper / paper and plastic / paper autoclave bags manufactured in accordance to EN 868. Loosely pack instruments in the bags to allow steam penetration to all instrument surfaces.

For STATIM 5000 G4:

The unwrapped instrument rack fitted with a maximum of 10 drying plates will hold 10 autoclave bags. Care must be taken to ensure that the combined weight of the loaded bags does not exceed 1.5 kg (3.3 lbs).

Rubber and Plastic Instruments

The following materials can be sterilized in the STATIM:

Nylon, polycarbonate (Lexan™), polypropylene, PTFE (Teflon™), acetal (Delrin™), polysulfone (Udel™), polyetherimide (Ultem™), silicone rubber, and polyester.

When loading rubber and plastic instruments in the tray, leave a space between the instruments and the cassette walls. This ensures that steam reaches all surfaces, and will promote drying.

4 Using Cassettes and Preparing Instruments

The following materials **cannot** be sterilized in the STATIM:

Polyethylene, ABS, styrene, cellulose, PVC, Acrylic (Plexiglas™), PPO (Noryl™), latex, neoprene, and similar materials.

Use of these materials may lead to instrument or equipment damage. If you are unsure of your instrument's material or construction, do not load into your STATIM until you have checked with the instrument manufacturer.

All Instruments

The STATIM is **NOT** intended for sterilizing textiles, liquids or biomedical waste. Instruments will remain sterile after a successful cycle until the cassette is disengaged from the unit. Unwrapped instruments, once exposed to ambient or external conditions, cannot be maintained in a sterile state. If sterile storage is desired, wrap the instruments to be sterilized in autoclave bags, according to the instrument manufacturer's instructions. Then allow the wrapped cycle to run until the air-dry phase is complete.

Best Practice Tip: Allow instruments (wrapped or unwrapped) to dry completely prior to handling. Wrapped or pouched instruments must not touch each other to promote drying and enable effective sterilization.

SciCan recommends users carefully choose the most appropriate sterilization cycle according to the recommendations of their leading infection control authorities and local regulatory guidelines / recommendations.

4.5 Using Biological and Chemical Indicators

Chemical process indicators suitable for steam sterilizers should be included in or on each package or load being sterilized. In addition, the weekly use of biological indicators, which allow you to ascertain whether the instruments have been exposed to sterilization conditions, is recommended.

4.6 Instrument Weight Guide

Instrument	Typical Instrument Weight
Scissors	30 g / 0.96 oz
Dental scalers	20 g / 0.64 oz
Forceps	15 g / 0.48 oz
Dental handpiece	40 to 60 g / 1.29 to 1.92 oz
Wrapped instrument rack	260 g / 8.35 oz
Unwrapped instrument rack	225 g / 7.23 oz
Suction cannula	10 g / 0.32 oz
Plastic mouth mirror	8 g / 0.25 oz
Impression tray	15 to 45 g / 0.48 to 1.45 oz
Plastic x-ray positioning ring	20 g / 0.64 oz

NOTE: The above weights are to be used as reference only. For exact weights of your instruments, consult the manufacturer's specifications.

5 Using your STATIM

Before using your STATIM for the first time, make sure the reservoir is full and the pump is properly primed. Refer to sections **3.4 Filling the Reservoir** and **3.5 Priming the Pump** for detailed instructions.

5.1 Selecting a Cycle

The STATIM 2000 G4 and 5000 G4 have seven sterilization cycles, each designed to sterilize using the parameters specified. Each cycle can be selected by pressing the UNWRAPPED, WRAPPED or RUBBER / PLASTIC cycle buttons.

The types of instruments, sterilization requirements and a graph depicting each cycle characteristics are described over the following pages. Load size requirements are listed in the 'Test Protocol' section.

5.1.1 Unwrapped Cycles

The STATIM 2000 G4 and 5000 G4 feature two 134°C type S and one 134°C type N UNWRAPPED sterilization cycles. At the end of the cycle's sterilization phase, air drying will commence for one hour. The type S cycles are used to sterilize solid and hollow metal instruments such as dental handpieces.


HOLLOW UNWRAPPED (S) 134°C / 3.5 min

HOLLOW UNWRAPPED (S) 134°C / 18 min

The type N cycle is used to sterilize solid metal instruments such as pliers, burrs, scalers and forceps.

SOLID UNWRAPPED (N) 134°C / 3.5 min

To select one of these cycles: Press the UNWRAPPED cycle button to scroll through the available cycles.

Once the desired cycle has been selected, press the **START** button.

The STATIM G4 will remember the last Unwrapped cycle type selected and will display it when you select the Unwrapped icon.


NOTE: If Process Enforced usage is enabled (where users must enter a PIN to Start and Stop a cycle), a PIN screen will appear after you have pressed START. Enter your PIN to start the cycle.

5 Using your STATIM


5.1.1 Unwrapped Cycles cont'd

HOLLOW UNWRAPPED (S) 134°C / 3.5 min

HOLLOW UNWRAPPED (S) 134°C / 18 min

SOLID UNWRAPPED (N) 134°C / 3.5 min

5 Using your STATIM

5.1.2 Hollow Wrapped Cycles


The STATIM 2000 G4 and 5000 G4 feature two WRAPPED 134°C type S sterilization cycles. These cycles are used to sterilize solid and hollow metal instruments which have been sealed in paper/paper or paper/plastic autoclave bags.

HOLLOW WRAPPED (S) 134°C / 3.5 min

HOLLOW WRAPPED (S) 134°C / 18 min

From the main menu, press the WRAPPED icon to scroll through the available cycles.


Once the desired cycle has been selected, press the **START** button.

The STATIM G4 will remember the last Hollow Wrapped cycle type selected and will display it when you select the Hollow Wrapped icon.

NOTE: A Process Challenge Device (PCD) is available for validation of HOLLOW WRAPPED (S) 134°C/3.5 min cycle.

HOLLOW WRAPPED (S) 134°C / 3.5 min

HOLLOW WRAPPED (S) 134°C / 18 min

5 Using your STATIM

5.1.3 Rubber and Plastics Cycles


The STATIM 2000 G4 and 5000 G4 feature two 121°C type S sterilization cycles used to sterilize solid unwrapped instruments constructed of metal or the materials listed in section 'Preparing and Loading Instruments'.

RUBBER / PLASTIC (S) 121°C / 15 min

RUBBER / PLASTIC (S) 121°C / 30 min

From the main menu, press the Rubber & Plastics icon to scroll through the available cycles.


Once the desired cycle has been selected, press the **START** button.

The STATIM G4 will remember the last Rubber / Plastics cycle type selected and will display it when you select the Rubber / Plastics icon.

RUBBER / PLASTIC (S) 121°C / 15 min

RUBBER / PLASTIC (S) 121°C / 30 min

5 Using your STATIM

5.1.4 Air Dry Only Cycle


This is not a sterilization cycle.

The load is considered sterile after the successful completion of the sterilization phase of the cycle. The Air Dry phase starts automatically after each sterilizing cycle and runs for 60 minutes.

Air drying may be interrupted by pressing the STOP button anytime after the sterilization phase of the cycle is complete. To ensure that the contents of the cassette are dry, the cycle should run for the full 60 minutes. Dryness is important for unwrapped instruments to prevent corrosion. For wrapped instruments, a dry wrap is required to maintain sterility.

If the STOP button is pressed during the air drying stage of the sterilization cycle, and the cassette has not been removed from the autoclave, the Air Dry Only cycle may be used to promote further drying. If the cassette has been removed from the autoclave, it may NOT be reinserted for the Air Dry Only cycle. If the cassette contains wrapped instruments and the wraps are not dry when the cassette is opened, the instruments must be handled in a manner that maintains their sterility and must be used immediately or resterilized.

NOTE: Sterilized instruments should only be handled once they are dry. Drying times can vary depending on the weight of the load, 60 minutes is based on a maximum load per cycle. If best practices are exercised (refer to section 'Preparing and Loading Instruments' and 'Maintenance') and the load is less than the maximum capacity, instruments may be dry in less than 60 minutes.

To start, press the Air Dry Only icon, then press the **START** button.

When started independently, this cycle will run for 60 minutes.

5 Using your STATIM

5.2 Running a Cycle

To operate each cycle, follow these steps.

1. Turn the power switch at the back of the unit to ON.

At start up, the unit will display the main menu.


2. Press the appropriate cycle button on the touchscreen to scroll through the available cycles.

The display will show the cycle name and parameters. The STATIM G4 features seven different cycles:

HOLLOW UNWRAPPED (S) 134°C / 3.5 min	or	HOLLOW WRAPPED (S) 134°C / 3.5 min	or	RUBBER / PLASTIC (S) 121°C / 15 min
HOLLOW UNWRAPPED (S) 134°C / 18 min		HOLLOW WRAPPED (S) 134°C / 18 min		RUBBER / PLASTIC (S) 121°C / 30 min
SOLID UNWRAPPED (N) 134°C / 3.5 min				

3. When you have found the cycle you require, press the START icon.

NOTE: If Process Enforced usage is enabled, a PIN screen will appear after you have pressed START. Enter your PIN to start the cycle.

5 Using your STATIM

5.2 Running a Cycle

When the cycle starts, the cycle parameters are displayed at the top of the screen.

Below it is the current phase. The unit's cycle counter is displayed at right.

A graph charts the progress of the cycle while current cycle information is displayed at right


As the cycle is running, various sounds will be heard. This is the normal functioning of the unit.

5 Using your STATIM

5.2 Running a Cycle cont'd


The buzzing noise during the air drying stage is the compressor operating. The air drying phase of the cycle may be interrupted at any time by pressing the STOP button.

When the automatic 60-minute air drying stage is finished and the sterilization cycle is successful, the touchscreen will display a Cycle Complete message and the reminder tone will sound until the STOP button is pressed or the cassette is removed from the unit.


Be careful. The metal parts will be hot and the cassette may contain hot steam.

5.3 Stopping a Cycle

To stop a cycle, press the **STOP** icon at the bottom right of the touchscreen. If the **STOP** button is pushed, the cassette is removed, or the unit detects a problem while operating, the cycle will stop. Once a cycle has been stopped, the **STOP** button must be pressed before another cycle can be started. The display will show:


If the display shows the message, **CYCLE FAULT** or **NOT STERILE**, the cassette contents are not sterile! See Section 9. Troubleshooting for more information.

If the air drying stage of the cycle is interrupted, do not store wrapped instruments that were in the cassette unless they are dry.

NOTE: If Process Enforced usage is enabled, a PIN screen will appear after you have pressed STOP. To release the load manually, press EN on the PIN screen. On the next screen, select **RELEASE LOAD MANUALLY**.

6 Storing and Retrieving Cycle Information

The STATIM G4 has an internal Datalogger capable of storing all cycle data on every cycle, whether successful or incomplete, for the lifetime of the unit. You can access this information through the touchscreen, through the web portal, using a USB storage device or by attaching a printer.

6.1 Retrieving Cycle Information Using the Touchscreen


1. From the main menu, press the USB icon.
2. The unit will record the last five successful cycles and the last five incomplete cycles. If you select a cycle from the list, it will display cycle information in a format similar to how it would be printed.
3. Use the arrow keys to scroll through and read.

NOTE: Regardless of whether you have a USB storage device attached to the unit or not, you can always see the last five successful cycles and the last five incomplete cycles. Use the STATIM G4 ConnectIM web portal to access all the cycle information stored on your STATIM from your computer. To connect your STATIM to a network consult the separate manual titled STATIM 2000/5000 G4 – ConnectIM User Guide.

6.2 Retrieving Cycle Information Using the USB Data Back Up

The USB storage device can be used to transfer cycle information stored in the unit to a computer. Best practice suggests this should be done once a week. To transfer data using the USB port, follow these steps:

1. Plug the USB storage device into the USB port.
2. The STATIM keeps track of what data has already been transferred to the USB storage device and will automatically load only new data.
3. When the activity light on the provided USB storage device stops blinking or the USB icon on the LCD turns from a flashing green to a solid grey, remove the USB storage device and transfer the information to your computer.

NOTE: If you select the USB storage device icon from the main menu, you will only be able to view the last five complete cycles and the last five incomplete cycles. To view all the cycles stored on the USB storage device, you must use your computer.

7 Printing Cycle Information

The STATIM G4 is equipped with an RS232 serial port to allow you to connect it to an external printer. (For a list of recommended printers, see the table below.)

Connecting to a Printer

To connect the printer, follow these steps:

1. Connect the external printer to the STATIM G4's RS232 port using the serial printer cable supplier with your printer.
2. Power on the printer.
3. On the STATIM G4, select  →  → 
4. Scroll to and select.
5. Use   to toggle to the Serial Printer option and select. Press  to save and move back to the Setup menu.

Adjusting your Print Settings

The STATIM G4 allows for several printer adjustments. You can access these settings from the user setup menu (see instructions above). Use the table below or your printer's operator manual to make the correct adjustments to your printer's , and

External Printers and Specifications

Suggested External Printers by SciCan	End Of Line CR/LF	Serial Port Bitrate	Printer user ° char
Epson TM-U220D(C31C515603)	CR/LF	9600	248 [0xF8]
Citizen IDP-3110-40 RF 120B	CR	9600	N/A
Star Micro SP212FD42-120	CR	9600	210 [0xd2]
Star Micro SP216FD41-120	CR/LF	9600	210 [0xd2]
Star Micro SP512MD42-R	CR/LF	9600	210 [0xd2]

8 Maintaining your STATIM

8.1 Cleaning the Cassette

Keeping the STATIM cassette clean is good clinical practice and assists in the proper functioning of the unit. SciCan recommends that the interior surface be cleaned at least once a week. Cleaning the inside of your cassette is very important if you regularly sterilize lubricated instruments.

1. Use dishwashing soap or a mild detergent that does not contain chlorine.
2. Scrub the inside of the cassette with a cleaning pad designed for use with Teflon™ coated surfaces.
3. After scouring, rinse thoroughly with water to remove all traces of the detergent.

Coating the entire inside surface with STAT-DRI Plus drying agent induces water to form an even coat on the inside surface, without beading. The water in contact with the hot cassette surfaces also evaporates much more efficiently. Spotting is minimized and instruments dry much better. STAT-DRI Plus should be applied every 10 cycles, and after every cassette cleaning.

STAT-DRI Plus is available from SciCan in a 2-ounce bottle (SciCan Part # 2OZPLUS) or an 8-ounce bottle (SciCan Part #8OZPLUS) or a 32-ounce bottle (SciCan Part #32OZPLUS).

8.2 Cleaning the Water Reservoir Filter

This water reservoir filter should be cleaned at least once a week or when required.

The filter can easily be removed and cleaned by placing the filter upside down under running water to wash away the particles until clean, and then placed back into the reservoir opening. If a replacement water reservoir filter is required, order part number 01-109300S.

8.3 Cleaning the Reservoir

Check the reservoir for dirt or particles. The reservoir may be cleaned by draining followed by cleaning and rinsing with steam process distilled water ONLY. Use of chemicals or cleaning agents is not recommended and could damage the unit.

8.4 Cleaning the Exterior Surfaces

Use a soft cloth moistened with soap and water to clean all exterior surfaces. Do not use harsh cleaning chemicals or disinfectants.

8.5 Changing the STATIM 2000 G4 Air Filter

The filter should be replaced every six months or after 500 cycles in order to maintain an adequate supply of clean air during the air drying cycle.

To change the filter, follow these steps:

1. Turn the power switch at the back of the unit **OFF**.
2. Remove and discard the old foam air filter **3**.
3. Install the new filter **3** (SciCan part no. 01-100207S).
4. Secure the filter plate **2** to the back of the compressor using the screw **1** retained during the disassembly procedure.


Figure 10



8 Maintaining your STATIM

8.6 Changing the Bacteria Retentive Air Filters

The filters should be replaced every six months or after 500 cycles to maintain an adequate supply of clean air during the air drying cycle.

To change the bacteria retentive air filter on the STATIM 2000 G4 and 5000 G4, follow these steps:

1. Power the STATIM **OFF**.
2. Disconnect tube A **1** from the bacteria retentive filter **2** and remove the filter from the filter bracket **3**. As you remove the filter from the bracket, note the orientation of the arrow mark on the filter.
3. When the filter is free of the bracket, carefully disconnect tube B **4** from the filter.
4. Before installing the replacement bacteria retentive filter **2** (SciCan order no. 01-102119S) check that the arrow mark on the filter matches the direction of the arrow on the bracket. Push the left hand filter fitting into tube B **4**.
5. Gently press the replacement filter into the filter bracket **3**. The arrow mark of the filter should be facing out and pointing to the left.
6. Re-connect tube A **1** to the right hand filter fitting.

8 Maintaining your STATIM


8.7 Replacing the Cassette Seal

To ensure optimum performance of your STATIM cassette autoclave, change the cassette seal every 500 cycles or every six months, whichever comes first. Replacement seals are available from SciCan (order number 01-100028S for the STATIM 2000 G4 and 01-101649S for the STATIM 5000 G4).

To change the cassette seal, follow these steps:

Place the cassette lid and the new seal on a clean work surface. Examine the position of the old seal in the cassette lid and arrange the new seal in the same orientation, next to the lid.

Remove the old seal and discard. Clean any residue out of the seal channel and flush out the channel with distilled water.

Lubricate the new seal with the liquid seal lubricant provided.

Insert the rounded edge of the seal under the round lip of the lid. Align the holes in the new seal with the holes in the lid.

NOTE: At every corner and at the holes in the lid, two square nibs should be visible. The nibs should fit flush with the lid's outer surface.

Ensure the seal is completely inserted. Feel around the periphery to ensure the seal is securely in place.

NOTE: During a cycle, steam may appear between the lid and the tray. If this persists, remove the cassette and check that the seal is correctly installed.

Be careful. The metal parts will be hot, and the cassette may contain hot steam.

Figure 13

8 Maintaining your STATIM

8.8 Maintaining Fluid Levels

1. The water reservoir level is continually monitored by your STATIM. If the reservoir is low, a red X will appear on the water icon  of the cycle select screen. Press the icon to go to the next screen to confirm that it is a water level issue and not a water quality issue.
2. If the reservoir is low, a  will appear next to WATER LEVEL.
3. To fill the reservoir, remove the cap from the top of the unit and fill the reservoir. We recommend using a funnel to minimize spills. Each time you refill the reservoir, empty the waste bottle and refill with water to the MIN line. Empty the waste bottle often to avoid unpleasant odors and discoloration of the contents. (A low-level chlorine-free disinfectant, prepared according to the manufacturer's instructions, may be added to the waste bottle to remedy this situation).

8.9 Reading Water Quality

1. The water quality is continually monitored by your STATIM. Use only steam-process distilled water containing less than 5 ppm total dissolved solids (having conductivity of less than 10 $\mu\text{S} / \text{cm}$) in the STATIM. If the water quality is beyond these tolerances, a red X will appear on the water icon  of the cycle select screen. Press the icon to go to the next screen to confirm that it is a water quality issue and not a water level issue.
2. If the water quality not suitable,  will appear next to micro S. and parts per million values.
3. Using the drain tube (see section 3.5 Priming the Pump, Figure 6) empty the contents of the reservoir into the water container and replace with steam-process distilled water containing less than 5 ppm total dissolved solids (having conductivity of less than 10 $\mu\text{S} / \text{cm}$).

8 Maintaining your STATIM

8.10 Preventative Maintenance Schedules

To ensure trouble-free performance, both the operator and the dealer must follow a preventative maintenance schedule.

NOTE: Please refer to your National, Regional, State or Safety laws for any additional reoccurring user testing that may be required.

The schedules below describe the necessary actions.

Operator		
Daily	Water Reservoir	<ul style="list-style-type: none"> • Drain after each working day.
	Waste Bottle	<ul style="list-style-type: none"> • Empty the waste bottle every time you refill the reservoir. Fill with tap water up to MIN line. • You can also add some chlorine-free disinfectant.
Weekly	Cassette	<ul style="list-style-type: none"> • Wash the interior of the cassette with chlorine-free dishwashing detergent or soap. Rinse thoroughly with water to remove all traces of the detergent and treat interior surfaces of the cassette with the STAT-DRI™ Plus drying agent to enhance the drying process. • Order more STAT-DRI™ Plus from SciCan quoting 2OZPLUS, 8OZPLUST, or 32OZPLUS.
	Water Filter	<ul style="list-style-type: none"> • Check the water reservoir filter every week and clean if necessary. Replace only if necessary.
	Biological and/or Air Filter	<ul style="list-style-type: none"> • Check the filter for dirt and moisture. Replace if dirty. Call for service if wet.
Every 6 months	Air Filter	<ul style="list-style-type: none"> • Replace every 500 cycles or six months (whichever comes first).
	Biological Air Filter	<ul style="list-style-type: none"> • Replace every 500 cycles or six months (whichever comes first).
	Cassette Seal	<ul style="list-style-type: none"> • Replace every 500 cycles or six months (whichever comes first).

Technician		
Once a year	Cassette	<ul style="list-style-type: none"> • Check the tray, lid and seal for damage. Replace if necessary.
	Biological Filter	<ul style="list-style-type: none"> • Inspect the biological filter for moisture.
	Solenoid Valve	<ul style="list-style-type: none"> • Inspect the valve and clean if dirty. Replace the plunger if defective.
	Pump	<ul style="list-style-type: none"> • Clean the filters, replace if dirty.
	Check Valve	<ul style="list-style-type: none"> • Remove the exhaust tube from the back of the unit during the air drying phase. Check for air coming from the fitting. • Remove the air compressor tube from the check valve inlet while running a cycle. Make sure no steam is leaking from the valve. Replace if there are any leaks.
	Water Reservoir	<ul style="list-style-type: none"> • Check the reservoir for dirt. Clean and rinse with steam process distilled water if necessary.
	Calibration	<ul style="list-style-type: none"> • Calibrate the unit.

8 Maintaining your STATIM

8.11 Shipping the Unit / Draining the Reservoir

Before you move the unit, you will need to drain the reservoir. To do so, follow these steps:

1. Place a water container below the unit.
2. Using the drain tube (see Section 3.5 Priming the Pump, Figure 6) empty the contents of the reservoir into the water container.
3. Remove any remaining water from the reservoir with a non-linting, absorbent towel.
4. Screw-in the three leveler feet found underneath the unit.
5. Repack the unit in the original packing materials and include all accessories originally included with the unit.
6. Specify heated and insured shipping.

9 Troubleshooting Your STATIM

Problem	Solution
<p>Unit does not power ON.</p>	<p>Check that the unit is plugged into a properly grounded outlet and that the power cord is firmly seated at the rear of the machine.</p> <p>Try another circuit. Power unit OFF for 10 seconds and then power ON again.</p> <p>Check the condition of the line circuit breaker or fuse.</p>
<p>There is water under the machine.</p> <div style="text-align: center;">   </div>	<p>Check that water was not spilled when refilling the reservoir. Make sure the plug in the exhaust tube is secured. Remove and reinsert the cassette. Attempt another cycle.</p> <p>Be careful. The metal parts will be hot, and the cassette will contain hot steam.</p> <p>The cassette is leaking. If water drips from the underside of the unit during operation, check the cassette seal for misalignment or damage and replace the seal if required.</p> <p>Be careful. The metal parts will be hot, and the cassette will contain hot steam.</p> <p>Attempt another cycle. If it still leaks attempt another cycle using a different cassette if possible.</p> <p>If the leak persists, turn the unit OFF, remove and unload the cassette, unplug the unit, and call your dealer.</p>
<p>Instruments do not dry.</p>	<p>Best drying occurs when the cycle continues to completion. Allow the cycle to finish. Make sure the instruments are loaded correctly in the cassette. Refer to section 4. Using Cassettes and Preparing Instruments.</p> <p>Check the unit leveling.</p> <p>Check air/biological filters and replace if dirty. Clean the inside of the cassette and treat with STAT-DRI Plus drying agent. Refer to section 8.1 Cleaning the Cassette. Examine the exhaust tube (tube to the waste bottle) for kinks.</p>

9 Troubleshooting con't.

Problem	Solution
<p>Instruments do not dry con't</p>	<p>If kinked, straighten the tube. If the tube cannot be straightened, remove it from the push-in fitting attached to the STATIM. Depress the collar on the fitting and with the other hand pull firmly on the tube. Once the tube is free of the fitting, cut the damaged section of tubing away using a sharp instrument. Be sure that you leave enough tube to reach the unit when you re-attach the tube to the exhaust fitting. If the tube is too short to remove a section, contact your SciCan dealer for a replacement.</p> <p>Make sure the compressor is working. To check, remove the exhaust tube from the waste bottle. Start the Air Drying Only Cycle, and place the free end into a glass of water. If there is not a strong, steady flow of bubbles, the compressor is not functioning properly. Contact your SciCan dealer.</p>
<p>Cycle interrupted — NOT STERILE, Cycle aborted — NOT STERILE and CYCLE FAULT messages.</p> 	<p>Wait a few minutes and attempt another cycle before proceeding to the next solution. Remove the cassette. Be careful. The metal parts will be hot and the cassette will contain hot steam. Inspect the cassette to ensure that the holes in the back of the seal are perfectly aligned, and that the flexible lip of the seal is completely free. Check the exhaust tube for kinks or obstructions. If kinked, straighten the tube.</p> <p>If the tube cannot be straightened, remove it from the push-in fitting attached to the STATIM. Depress the collar on the fitting and, with the other hand pull firmly on the tube. Once the tube is free of the fitting, cut the damaged section of tubing away using a sharp instrument. Be sure that you leave enough tube to reach the unit when you re-attach the tube to the exhaust fitting. If the tube is too short to remove a section, contact your SciCan dealer for a replacement.</p> <p>Check that the STATIM has not inadvertently been exposed to any electrical interference. Refer to section 3.1 Positioning your Unit.</p> <p>Try running another cycle. If the problem persists, record the cycle fault message number and contact your dealer.</p>

9 Troubleshooting con't.

Problem	Solution
<p>Excessive steam issuing from the front of the machine.</p> 	<p>Remove and reinsert the cassette. Attempt another cycle. Remove and check the cassette seal for misalignment or damage. Replace the seal if required. Be careful as the metal parts will be hot and the cassette will contain hot steam.</p> <p>If the leak persists, turn the unit OFF, remove and unload the cassette and contact your SciCan dealer.</p>
<p>Machine will not start and touchscreen shows:</p> 	<p>Press on the icon to confirm whether it is a water level problem or a water quality problem. If it is a water quality problem, you have likely used water which is not steam-process distilled or is improperly distilled.</p> <p>Empty the reservoir and refill with steam-process distilled water containing less than 5 ppm total dissolved solids (having conductivity of less than 10 $\mu\text{S} / \text{cm}$). If you have the water conductivity meter, check the quality of the water before refilling the reservoir. To empty the reservoir, see section 8.11 Shipping the unit / Draining the Reservoir.</p>
<p>Machine will not start and touchscreen shows:</p> 	<p>Press on the icon to confirm whether it is a water level problem or a water quality problem.</p> <p>If the level of the water in the reservoir is low, refill the reservoir. Refer to the steps described in section 3.4 Filling the Reservoir.</p>
<p>The printer does not work.</p>	<p>Make sure that the printer cable is connected securely with the connector on the back of the STATIM. Make sure that the printer is powered ON. Power unit OFF for 10 seconds and then power ON again.</p>
<p>Time and date are incorrect.</p>	<p>The time and date have not been set. See Section 3. Setting up your STATIM, for time and date instructions.</p>

9 Troubleshooting con't.

Problem	Solution
<p>Touchscreen is blank/white</p>	<p>Power was interrupted during a firmware upgrade. Power off the unit and power it on again. It can take up to 6 minutes before the main menu screen appears.</p>
<p>Touchscreen is blank/ dark</p>	<p>Check power source.</p>
<p>USB storage device does not contain the last print out</p>	<p>Re-insert the USB storage device and wait for the data to copy over again. If problem persists, back up all the information you have on the USB device and reformat it. NOTE: You can always access all your unit's cycle information through the unit's web portal.</p>
<p>Touchscreen shows:</p> 	<p>An X over the connectivity icon means the unit is not connected to a network. If it is supposed to be connected to a network and the X is visible, it is because the unit is unable to acquire an IP address. To resolve the issue, try some of the following:</p> <ul style="list-style-type: none"> ● Check that the router is functioning properly ● Check the LAN cable (try a new cable if possible) ● Make sure your router assigns IP addresses automatically. ● Renew the IP address by following these steps: <ol style="list-style-type: none"> 1. Scroll through the setup menu to NETWORK SETUP and select. 2. Select RENEW IP
<p>Unit is not sending emails</p>	<p>Check email settings by using the TEST button on the unit's web portal. From the SETUP web page, select the TOOLS tab. Click on TEST to check your router, unit, and Internet connections. If all settings appear to be OK. Go to the unit's touchscreen and renew the IP address by following these steps:</p> <ol style="list-style-type: none"> 1. Scroll through the setup menu to NETWORK SETUP and select. 2. Select RENEW IP
<p>Not receiving emails from the unit</p>	<p>Check your spam filter. Be certain the unit has been identified as an accepted email source. Ensure that you have accepted the SciCan Privacy policy by putting a check mark in the box on the CONTACTS page of your web portal.</p>

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1.7 mm dia x 104 mm length #10-0008-00	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Trocac sleeve, arthroscopy, 2 rotating stop cocks,	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Trocac sleeve, hysteroscopy diagnostic sheath, 1 fixed stop cock, 2.7 mm dia x 302 mm length #10-0049-00	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min

Alcon ophthalmic handpiece

NeoSonix Phaco handpiece	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
--------------------------	--

Miltex medical instruments

Frazier needle 26-778	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Yeoman biopsy forceps with rotating shaft 28-304	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Kerrison rongeur 18-1994	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Frazier-Ferguson tube 19-570	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Yankauer suction tube 2-104SS	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Menghini biopsy needle 13-150	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min

Becton Dickinson

Needle, 30G1	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
--------------	--

Medical Workshop

Hockeystick Forceps, membrane peeling mw-1925	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
---	--

- * 1. The PCD (Process Challenge Device) was designed exclusively for use in STATIMs that are compliant with EN13060.
2. The PCD is not for use with STATIMs which are not in compliance with EN13060 or with the extended cassette.
3. PCD specification available upon request.
4. Hollow load A test applicable – rationale available upon request.

10 Test Protocol con't.

Cycle	STATIM 2000 G4	STATIM 5000 G4
	Cold Unit With Max. Load + Drying Phase	Cold Unit With Max. Load + Drying Phase
	Warm Unit Without Load + Drying Phase	Warm Unit Without Load + Drying Phase
SOLID UNWRAPPED (N) 134°C / 3.5 min	9:15 + 60:00	13:15 + 60:00
	6:45 + 60:00	8:45 + 60:00
HOLLOW UNWRAPPED (S) 134°C / 3.5 min	11:45 + 60:00	17:30 + 60:00
	8:05 + 60:00	10:50 + 60:00
HOLLOW UNWRAPPED (S) 134°C / 18 min	26:15 + 60:00	32:00 + 60:00
	22:35 + 60:00	25:20 + 60:00
HOLLOW WRAPPED (S) 134°C / 3.5 min	15:35 + 60:00	24:00 + 60:00
	10:40 + 60:00	15:30 + 60:00
HOLLOW WRAPPED (S) 134°C / 18 min	30:05 + 60:00	38:30 + 60:00
	25:10 + 60:00	30:00 + 60:00
RUBBER / PLASTIC (S) 121°C / 15 min	20:15 + 60:00	22:50 + 60:00
	18:40 + 60:00	20:20 + 60:00
RUBBER / PLASTIC (S) 121°C / 30 min	35:15 + 60:00	37:50 + 60:00
	33:40 + 60:00	35:20 + 60:00

11 Ordering Spare Parts

STATIM G4 SPARES	
01-100028S	Cassette Seal (2000)
01-101649S	Cassette Seal (5000)
01-103865S	Seal Lubricant
01-100207S	Compressor Filter (2000)
01-102119S	Filter Biological
01-101783S	Reservoir Cap and Filter
01-109300S	Water Reservoir Filter Kit
01-100204S	Exhaust Tube
01-100724S	Condenser Bottle w / o Condenser
01-100735S	Waste Water Bottle Fitting
01-100780S	Bumper
01-100812S	Condenser Bottle
01-104093S	Exhaust Tube 3 m long
01-104343S	Plug - Drain Tubing
01-108340S	Statim PCD – spare parts
01-100782S	Push-In Fitting (2000)
01-101755S	Push-In Fitting (5000)
01-103945S	Rack-Tray Unwrapped Instr. Kit (STATIM 2000)
01-112409S	Cassette Lid (2000 G4)
01-112410S	Cassette Handle - Lid (2000 G4)
01-112386S	Cassette Lid (5000 G4)
01-112387S	Cassette Handles Tray/Lid (5000 G4)
01-112388S	Cassette Handle - Lid (5000 G4)
01-112511S	Cassette Lid (5000 Ext G4)
01-112512S	Cassette Handles Tray/Lid (5000 Ext G4)
01-112513S	Cassette Handle - Lid (5000 Ext G4)
01-103557S	Power Cord Replacement Denmark (det.)
01-101766S	Power Cord UK
01-101768S	Power Cord Switzerland
01-101769S	Power Cord Italy
01-101779S	Power Cord Europe

STATIM G4 ACCESSORIES	
01-101709S	Mesh Rack (5000)
01-106653	Mesh Rack - STATIM 2000
01-112408S	Cassette Tray (2000 G4)
01-112407S	Cassette Tray with mesh rack (2000 G4)
01-112406S	Cassette Complete (2000 G4)
01-112510S	Cassette Tray (5000 Ext G4)
01-112509S	Cassette Complete (5000 Ext G4)
01-112385S	Cassette Tray (5000 G4)
01-112384S	Cassette Complete (5000)
01-106325	Container Endoscope Complete (STATIM 5000)
01-103935	STAT-DRI Plates (Qty 5) STATIM 5000
01-103923	Condenser Additional Bottle
2OZPLUS	STAT-DRI Plus 2 oz.
8OZPLUST	STAT-DRI Plus 8 oz.
32OZPLUS	STAT-DRI Plus 32 oz.
99-108332	Chemical Emulators (Class 6) 134°C/3.5min)
01-108341	PCD Final Assembly Kit

12 Warranty

Limited Warranty

For a period of one year, **SciCan** guarantees that the **STATIM 2000 / 5000 G4**, when manufactured by **SciCan** in new and unused condition, will not fail during normal service due to defects in material and workmanship that are not due to apparent abuse, misuse, or accident.

The one year warranty will cover the performance of all components of the unit except consumables such as the cassette seal, the compressor filter and the microbiological filter, provided that the product is being used and maintained according to the description in the user's manual.

A two-year warranty will be applied specifically to the water pump, steam generator, and printed circuit board (PCB) provided that the product is being used and maintained according to the description in the user's manual.

In the event of failure due to such defects during this period of time, the exclusive remedies shall be repair or replacement, at **SciCan's** option and without charge, of any defected part(s) (except gasket), provided **SciCan** is notified in writing within thirty(30) days of the date of such a failure and further provided that the defective part(s) are returned to **SciCan** prepaid.

This warranty shall be considered to be validated, if the product is accompanied by the original purchase invoice from the authorized **SciCan** dealer, and such invoice identifies the item by serial number and clearly states the date of purchase. No other validation is acceptable.

After one year, all **SciCan's** warranties and other duties with respect to the quality of the product shall be conclusively presumed to have been satisfied, all liability therefore shall terminate, and no action or breach of any such warranty or duty may thereafter be commenced against **SciCan**.

Any express warranty not provided hereon and any implied warranty or representation as to performance, and any remedy for breach of contract which, but for this provision, might arise by implication, operation of law, custom of trade or course of dealing, including any implied warranty of merchantability or of fitness for particular purpose with respect to all and any products manufactured by **SciCan** is excluded and disclaimed by **SciCan**.

If you would like to learn more about **SciCan** products and features or register your warranty online, visit our website at www.scican.com.

13.1 STATIM 2000 G4

Machine Dimensions:	Length:	49.5 cm
	Width:	41.5 cm
	Height:	15 cm
Cassette Size (External):	Length:	41 cm (includes handles)
	Width:	19.5 cm
	Height:	4 cm
Cassette Size (Internal):	Length:	28 cm
	Width:	18 cm
	Height:	3.5 cm
Sterilization Chamber Volume:		1.8 L
Reservoir Volume:		4.0 L
Weight (Without water):		22 kg
Clearance required:	Top:	5 cm
	Sides:	5 cm
	Back:	5 cm
	Front:	48 cm
Minimum charge in the water reservoir:		550 mL
PRV value (pressure relief valve):		43.5 PSI
Electrical Rating:		220 - 240 V, 50 Hz, 6 A
Ethernet Port:		10/100 Base-T
USB Port:		USB 2.0
Current:		AC
Protection Class:		I
Protection:		covered
Ambient Operating Temperature:		5°C - 40°C
Sound levels:		Mean - 56 dB, Peak - 65 dB
Humidity:		80% Max.
Max. Altitude:		2000 m

13.2 STATIM 5000 G4

Machine Dimensions:	Length:	60 cm
	Width:	41 cm
	Height:	19 cm
Cassette Size (External):	Length:	49.5 cm (includes handles)
	Width:	19.5 cm
	Height:	8 cm
Cassette Size (Internal):	Length:	38 cm
	Width:	18 cm
	Height:	7.5 cm
Sterilization Chamber Volume:		5.1 L
Reservoir Volume:		4.0 L
Weight (Without water):		34 kg
Clearance required:	Top:	5 cm
	Sides:	5 cm
	Back:	5 cm
	Front:	57 cm
Minimum charge in the water reservoir:		550 mL
PRV value (pressure relief valve):		43.5 PSI
Electrical Rating:		220 - 240 V, 50 Hz, 6 A
Ethernet Port:		10/100 Base-T
USB Port:		USB 2.0
Current:		AC
Protection Class:		I
Protection:		covered
Ambient Operating Temperature:		5°C - 40°C
Sound levels:		Mean - 57 dB, Peak - 65 dB
Humidity:		80 % Max.
Max. Altitude:		2000 m

STATIM 2000/5000 G4

CASSETTE AUTOCLAVE™

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SciCan | Your Infection Control Specialist™

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Make - Model	Cycle
Dental instruments	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Kavo - GENTLE force 7000C	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Kavo - Gentle PowerLux 2.5 LPA	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
W & H - Trend LS, W D-56	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
W & H - Trend HS, TC-9SRM	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Kavo Super-Torque LU X/640 B	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Kavo - NTRAMatic LU X3, 20 LH	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
NSK - FANA AX	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
NSK - ATL118040	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
NSK - TM AX	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
STAR-430 SW L	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Sirona - T1 Classic, S 40 L	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Sirona - T1 Control, TC3	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
M - Westfaldhörn	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Beh - Ax - Boel L	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
W & H - WS-75	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
W & H - WA-99 LT	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
W & H - TA-98 LC	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
B & L Ophthalmology instruments	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Globe irrigating cannula 30g E 4894	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Lask cannula E 4989	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Gills irrigating-aspating cannula E 4932	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Niham in hydrosection cannula 26g E 4421 H	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Irrigating-aspating handpiece MVS 1063C	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Rudolf Medicaltechnik GmbH Endoscope accessories	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Trocars beve, arthroscopy, 2 rotating stop cocks, 1.7 mm dia x 104 mm length #10-4008-00	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Trocars beve, hysteroscopy diagnostic sheath, 1 fixed stop cock, 2.7 mm dia x 302 mm length #10-0049-00	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Alcon ophthalmic handpiece	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Neosonic Phaco handpiece	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Milux medical instruments	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Frazier needs 26-778	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Yeoman biopsy forceps with rotating shaft 28-304	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Kerison rongeur 18-1994	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Frazier Ferguson tube 19-570	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Yankauer suction tube 2-104SS	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Menghini biopsy needs 13-150	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Becton Dickinson	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Needle, 30G1	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Medical Workshop	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min
Hockeystick Forceps, membrane peeling m w-1925	HOLLOW / UNWRAPPED (S) 134°C / 3.5 min

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STATIM G4 SPARES	
01-100028S	Cassette Seal (2000)
01-101649S	Cassette Seal (5000)
01-103865S	Seal Lubricant
01-100207S	Compressor Filter (2000)
01-102119S	Filter Biological
01-101783S	Reservoir Cap and Filter
01-109300S	Water Reservoir Filter Kit
01-100204S	Exhaust Tube
01-100724S	Condenser Bottle w / o Condenser
01-100735S	Waste Water Bottle Fitting
01-100780S	Bumper
01-100812S	Condenser Bottle
01-104093S	Exhaust Tube 3 m long
01-104343S	Plug - Drain Tubing
01-108340S	STATIM PCD – spare parts
01-100782S	Push-In Fitting (2000)
01-101755S	Push-In Fitting (5000)
01-103945S	Rack-Tray Unwrapped Instr. Kit (STATIM 2000)
01-112409S	Cassette Lid (2000 G4)
01-112410S	Cassette Handle - Lid (2000 G4)
01-112386S	Cassette Lid (5000 G4)
01-112387S	Cassette Handles Tray/Lid (5000 G4)
01-112388S	Cassette Handle - Lid (5000 G4)
01-112511S	Cassette Lid (5000 Ext G4)
01-112512S	Cassette Handles Tray/Lid (5000 Ext G4)
01-112513S	Cassette Handle - Lid (5000 Ext G4)
01-103557S	Power Cord Replacement Denmark (det.)
01-101766S	Power Cord UK
01-101768S	Power Cord Switzerland
01-101769S	Power Cord Italy
01-101779S	Power Cord Europe

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Rudolf Medizintechnik GmbH Endoscopy accessories