

velopex[®] AquaCare
INTERNATIONAL

TOTAL CONTROL OF YOUR

POWDER POWER

MINIMUM EFFORT

MAXIMUM CLEAN


 **optident[™]**

 **HENRY SCHEIN[®]**
DENTAL

ACCREDITED BY
STYLEITALIANO[™]
SHADE GUIDES

EDUCATED & ACCREDITED BY STYLEITALIANO

STEP 4 IN THE DIRECT WORKFLOW



"CLEAN TOOTH TISSUE IS ESSENTIAL TO BOND TO FOR THE MOST PREDICTABLE DIRECT AND INDIRECT RESTORATIONS; THANKS TO THE EFFICIENCY AND BRILLIANCE OF AQUACARE, AN EVERYDAY TOOL THAT IS HARD TO LIVE WITHOUT IN MY PRACTICE."

AJAY JUNEJA - STYLEITALIANO HONORARY MEMBER



HOW CLEAN IS CLEAN?

THE UNIT OF CHOICE BY WORLD LEADING CLINICIANS

INNOVATION MAKES IT UNLIMITED

RESEARCH AND DEVELOPMENT IS OUR LIFEBLOOD. OUR DENTAL UNDERSTANDING, APPRECIATION AND EXCELLENCE IS TRANSFORMING THE PARTICLE ABRASION FIELDS OF RESTORATIVE, PROPHYLAXIS AND CLEAN DENTISTRY.

TESTED, APPROVED & USED DAILY
BY KEY OPINION LEADERS AROUND THE WORLD



AQUACARE ALLOWS ME FINALLY TO SANDBLAST OR PROPHYLAXIS INSIDE AND OUTSIDE THE MOUTH, WITH OR WITHOUT WATER, WITH THE EXACT AMOUNT OF MEDIA I WANT TO DELIVER. I'M IN FULL CONTROL NOW OF THOSE IMPORTANT PROCEDURES."

DR. WALTER DEVOTO, ITALY, STYLEITALIANO



CLICK TO WATCH **KATHERINE LOSADA**
PREPARE AND CLEAN A RESTORATION
WITH AQUACARE

WITH MORE POWDERS IN DEVELOPMENT AND THE ABILITY TO CHANGE IN JUST ONE SECOND, THE AQUACARE READILY ALLOWS YOU TO EMBRACE ABRASION, CLEANING, DESENSITISING, REMINERALISING...

EXPLORE THE VERSATILE CAPABILITIES



PERIO



ORTHO




ENDO



PEDIATRIC



IMPLANTS



HYGIENE



ADHESIVE DENTISTRY / BONDING / CLEANING / DESENSITISING
LABORATORY / LITHIUM DISILICATE / MI / PEKKTON®
REMINERALISING / RESTORATIVE DENTISTRY
STAIN REMOVAL / VENEERS / ZIRCONIA OXIDE...
SUPRA OR SUB GINGIVALLY

ACCREDITED BY
STYLEITALIANO
SHADE GUIDES


PATENTED DELIVERY SYSTEM...

...VIA VENTURI SINGLE-USE TIPS, OUR AQUASOL SHROUD IS MOVING PARTICLE ABRASION FROM ITS DUSTY PAST TO ITS GLORIOUS FUTURE.

VIRTUALLY UNLIMITED, AQUACARE WITH ITS UNIQUE CARTRIDGE ALLOWS SIMPLE SWITCHING TO SUIT THE CLINICIAN'S TREATMENTS. IS AQUACARE FUTURE PROOF?

NEW POWDERS, NO PROBLEM, NEW SOLUTIONS, NO PROBLEM.

CASE STUDIES



"CLEANING AND CONDITIONING OLD COMPOSITES WITH AQUACARE 50 MICRON, AND INSERTION OF 2 UNICA MATRICES IN THE CENTRALS."

DR JORDI MANAUTA, ITALY, STYLEITALIANO



"A BIT OF SYLC ON LOW PRESSURE CLEANS UP STAINING NICELY AROUND BRACKETS."

DR THOMAS TAHA, UNITED KINGDOM


FURTHER CASE STUDIES



"AIR POLISHING WITH A SODIUM BICARBONATE TO REMOVE STAINING."
PROF. PEET VAN DER VYVER, SOUTH AFRICA



*"AIR ABRASION FROM BUCCAL AND LINGUAL TO REMOVE ALL DEBRIS.
I USE AN AQUACARE."*
DR LINCOLN HARRIS, AUSTRALIA



*"MINIMALLY INVASIVE CAVITY PREPARATION OF
A CARIOUS LESION ON THE OCCLUSAL SURFACE A
MANDIBULAR THIRD MOLAR USING AQUACARE 29µM
ALUMINIUM OXIDE AIR ABRASION POWDER."*
DR LOUIS MACKENZIE, UNITED KINGDOM



*"PARTICLE ABRASION WITH 53 MICRON ALUMINIUM OXIDE TO
REMOVE BONDING RESIN USED FOR TEMPORARIES, TO CLEAN
ADHESIVE SURFACES AND IMPROVE BOND STRENGTHS BEFORE
CEMENTING FOUR DISILICATE VENEERS."*
DR LORENZO VANINI, ITALY



*"PREPARATIONS AFTER COMPOSITE RESIN BUILD-UPS TREATED
FOR ADHESIVE CEMENTATION. CLEANING HAS BEEN MADE
USING AQUACARE DEVICE. POWDER PROCUT 29 MICRONS
AL₂O₃ WITH LOW PRESSURE.."*
PROF DR FRANCESCO MANGANI, ITALY

- AQUACARE INCREASES BOND STRENGTH WITH ITS TEXTURING OF THE CAVITY SURFACE
- FASTER, MORE EFFICIENT CLEANING FOR REMOVING STAINS FROM ALL TOOTH ASPECTS, EVEN ACCESSIBLY DIFFICULT PITS AND CREVICES,WHILST ALSO TREATING SENSITIVITY.
- MINIMAL NEED FOR ANAESTHETICS COMBINED WITH THE LACK OF VIBRATION, HEAT GENERATION AND NOISE LEADS TO REDUCED PATIENT ANXIETY.
- INCREDIBLY VERSATILE ACROSS A WIDE VARIETY OF RESTORATIVE DENTAL PROCEDURES.
- CONVENIENT, EASY TO USE COLOUR-CODED CARTRIDGE SYSTEM MAKES SWAPPING BETWEEN TYPES OF TREATMENT EFFORTLESS.

AQUACARE STAND

THE AQUACARE STAND HAS BEEN DESIGNED TO BE COMPACT, SECURE AND MOBILE. THE SUPPORT SHELF ENSURES THE UNIT CANNOT SLIP WHEN THE STAND IS MOVED.

THE STAND ENABLES THE AQUACARE UNIT TO BE MOVED CONVENIENTLY INTO PLACEWHEN REQUIRED.

THE STAND IS AVAILABLE WITH OR WITHOUT THE HIGHLY EFFICIENT AIR DRYER

STAND (NON-DRYAIR) CODE: 1196090

STAND (DRYAIR) CODE: 1196091



WHY OUR UNITS & STANDS?

- INNOVATIVE
- VERSATILE
- UNIQUE
- UNLIMITED

AQUACARE DENTAL AIR ABRASION UNIT

AQUACARE IS A DENTAL AIR ABRASION UNIT, ENABLING YOU TO PROVIDE CONTACT-FREE AND MINIMALLY INVASIVE DENTISTRY (MID), MAKING ROUTINE DENTAL PROCEDURES PAINLESS FOR PATIENTS.

AQUACARE USES BREAKTHROUGH TECHNOLOGY TO DELIVER A COMFORTABLE, QUICK CLEAN PROCEDURE, BY USING A FINE STREAM OF FLUID COMBINED WITH A TINY VOLUME OF POWDER DIRECTED AT THE TEETH TO BE TREATED.

HEIGHT: 185MM / WIDTH: 265MM / DEPTH: 240MM / WEIGHT: 3.066KG



ALSO AVAILABLE - AQUACARE SINGLE UNIT

AIR ABRASION UNIT CODE: 1167301

SINGLE UNIT CODE: 1167300

WHY OUR POWDERS?

- SUPERDRY POWDER
- BUILT-IN FLOW CONTROL
- QUALITY ASSURED
- MULTIDISCIPLINARY RANGE
- USER-FRIENDLY WITH EASY COLOUR CODE CARTRIDGES


WE'VE DESIGNED AQUACARE SO THAT IT WORKS IDEALLY AND ALLOWS FOR THE BROADEST RANGE OF DENTAL CASES. THAT'S WHY WE INVENTED A CARTRIDGE SYSTEM: THEY'RE EASY TO SWITCH, REPLACE AND AS MORE POWDERS BECOME AVAILABLE THEY ARE FUTURE PROOF.



29 µm ALUMINIUM OXIDE

AL₂O₃ ABRASIVE POWDER, FOR ROUTINE DENTAL PROCEDURES.

4PK CODE: 1083618



53 µm ALUMINIUM OXIDE

AL₂O₃ ABRASIVE POWDER, FOR HYPOMINERALISE AND HEAVIER DENTAL PROCEDURES.

4PK CODE: 1083617



SODIUM BICARBONATE SOFT

MINT FLAVOUR AND WITH OVER 40% REDUCTION IN PARTICLE SIZE FOR MAXIMUM PATIENT COMFORT.

4PK CODE: 1197356



SODIUM BICARBONATE

NAHCO₃ CLEANING POWDER, EFFECTIVE IN THE REMOVAL OF THE BIOFILM AND GENERAL PROPHY CLEANING.

4PK CODE: 1109968



SYLC DESENSITISING

BAG45S5 BIOGLASS CLEANING POWDER REMINERALISING AND DESENSITISING, STAIN REMOVAL FOR SUPRA & SUB GINGIVAL.

4PK CODE: 1161872

WHY OUR RANGE OF TIPS?

■ PRECISION ■ EXCELLENCE ■ INNOVATION ■ DETAIL-ORIENTED

AQUACARE PROVIDE 3 TYPES OF TIPS TO HELP YOU WITH YOUR CHOICE OF PROCEDURES AND TREATMENT.

EACH OF OUR TIPS ADAPT TO PROCEDURES AND TREATMENTS, HELPING YOU TO REACH 99% PRECISION.



AQUACARE STANDARD TIP:

SUITABLE FOR ALL POWDERS AND BOTH HANDPIECES.

50PK CODE: 1083616

AQUACARE iTIP:

DESIGNED FOR INTERPROXIMAL AND DIASTEMA USE, THOUGH ROUTINELY USED IN MANY OTHER APPLICATIONS ONLY TO BE USED WITH THE 0.6MM HANDPIECE.

50PK CODE: 1191840

AQUACARE PLUS TIP:

IDEAL FOR LARGER PARTICLE SIZE AND WIDER COVERAGE, E.G. HYPOMINERALISED CASES & HYGIENE USE.

50PK CODE: 1191839



0.6MM
SILVER
HANDPIECE:

CODE: 1191845

■ LIGHT WEIGHT
■ PRECISE CONTROL



0.8MM
GOLD
HANDPIECE:

CODE: 1191846

■ AUTO-CLAVE
■ NEVER CLOG

DISCOVER AQUACARE HANDPIECES

THE PATENTED HANDPIECES DELIVER A HIGH-SPEED FLUID STREAM, SUPPLIED FROM AN INVERTED BOTTLE DELIVERY SYSTEM, THE AQUASOL, TO CREATE A LIQUID CURTAIN AROUND THE ABRASIVE MEDIUM.

WHY AQUACARE AQUASOL?

- KEEPS THE MESS AWAY
- VANILLA TASTE FOR PATIENT COMFORT
- DESIGNED WITH A FINE AMOUNT OF ETHANOL (17.5%) TO AUTO-CLEAN YOUR AQUACARE WHILE WORKING (PREVENTS THE BUILD-UP OF BACTERIA)
- SUITABLE FOR ALL AQUACARE UNITS



AQUASOL FLUID IS IDEAL FOR USE WITH THE AQUACARE AS THE FLUID COMPONENT TO AQUAABRASION™. WITH A BEAUTIFUL VANILLA FRAGRANCE, YOUR PATIENTS WILL ENJOY THEIR AQUACARE TREATMENT. A 17.5% ETHANOL CONTENT HELPS TO STERILISE THE INNER WORKINGS OF YOUR AQUACARE, HELPING TO MAINTAIN HYGIENIC PROCEDURES.





















1x 500ML BOTTLE CODE: 1109511

6x 500ML BOTTLES CODE: 1191841

AQUACARE PRICING

1X AIRLINE IS INCLUDED IN INSTALLATION. FOR ADDITIONAL AIRLINES CHARGES APPLY.

	AQUACARE - JDE 1167301 RRP £3245.00		AQUACARE HANDPIECE SILVER 0.6MM JDE 1191845 RRP £122.97		AQUACARE SODIUM BICARBONATE SOFT POWDER JDE 1197356 RRP £57.87
	AQUACARE SINGLE - JDE 1167300 RRP £2165.00		AQUACARE HANDPIECE GOLD 0.8MM JDE 1191846 RRP £122.97		AQUASOL FLUID - JDE 1109511 RRP £14.47
	AQUACARE STAND - JDE 1196090 RRP £245.94 AQUACARE STAND (DRY AIR VERSION) - JDE 1196091 RRP £475.00		AQUACARE 29 MICRON ALUMINIUM OXIDE POWDER JDE 1083618 RRP £64.53		AQUASOL FLUID ECONOMY PACK - JDE 1191841 RRP £66.31
	AQUACARE STANDARD TIPS - JDE 1083616 RRP £24.74		AQUACARE 53 MICRON ALUMINIUM OXIDE POWDER JDE 1083617 RRP £64.53		AQUACARE HANDLES - JDE 1191847 RRP £19.29
	AQUACARE PLUS TIPS - JDE 1191839 RRP £24.74		AQUACARE SODIUM BICARBONATE JDE 1109968 RRP £53.05		AQUACARE DRYAIR UNIT - JDE 1191843 RRP £229.06
	AQUACARE ITIPS - JDE 1191840 RRP £24.74		AQUACARE SYLC POWDER JDE 1161872 RRP £129.63		DENTAL AIR ABRASION CABINET - JDE 1191842 RRP £245.94

FOR MORE INFORMATION CLICK BELOW:

CALL THE SALES TEAM ON
+44 (0) 1943 60 50 50

VISIT OUR WEBSITE
AND USE OUR LIVE CHAT

BOOK A FREE
ONLINE CONSULTATION

REQUEST A
CALL BACK

 **optident**[™]

 **HENRY SCHEIN**[®]
DENTAL

OPTIDENT.CO.UK / SALES@OPTIDENT.CO.UK /

 [@OPTIDENTLTD](https://www.instagram.com/optidentltd) /  [OPTIDENTLTD](https://www.facebook.com/optidentltd)

**BECAUSE
INGREDIENTS
MATTER**